

# Guide to Time in Range (TIR) and Coefficient of Variation (CV) Visualizations

## Daily Views

### Time in Range / Coefficient of Variation

Daily percentages for Time in Range and CV show what happened on each specific day.

- These values don't explain *why* they occurred.
- For example, 70 percent in range with a higher CV could reflect afternoon spikes, overnight stability, or a mix of both.
- Daily summaries highlight day-to-day differences.
- Only multi-day, hour-by-hour views reveal the underlying patterns that create those results.

## Hourly Views by Day

### Time in Range by Hour

Hourly Time in Range provides more context than a single daily value.

- On 11/28, thirteen hours stayed within 70 to 140 mg/dL.
- Hours after lunch and dinner lowered the overall percentage.
- This view identifies when glucose moved out of range.
- The pattern still varies day by day, and single-day views do not show the consistent trends that only appear when many days are combined.

### Coefficient of Variation by Hour

Hourly CV shows how much glucose fluctuated throughout the day.

- On 11/28, variability stayed low overnight and through the morning.
- Variability increased in the afternoon and early evening.
- This highlights when glucose became more variable, but the pattern still changes from one day to the next.
- A single day's hourly CV provides helpful context, yet repeatable trends only emerge when many days are viewed together.

## Multi-Day Hourly Views

### Time in Range by Hour - Multi-Day Views

Combining 7, 14, 30, or 90 days smooths out day-to-day variation.

- This reveals the underlying structure of the day rather than the noise of any single day.
- Individual days may differ, but a consistent pattern emerges.
- The typical rhythm in this case shows steady overnight hours, mid-day dips, and a return to range in the evening.

### CV by Hour - Multi-Day Views

Multi-day CV reduces noise present in single-day readings.

- It highlights the underlying pattern of glucose variability.
- Individual days may show sharp rises or dips.
- When 7, 14, 30, or 90 days are combined, a consistent trend becomes clear: low variability overnight, rising variability through the afternoon, and settling in the evening.
- This view shows the typical daily variability pattern rather than fluctuations from any single day.